OVITEX® CORE

REINFORCED TISSUE MATRIX

OVITEX® 1S

REINFORCED TISSUE MATRIX

OVITEX® 2S

REINFORCED TISSUE MATRIX

POLYMER REINFORCEMENT

PERMANENT
Polypropylene

SCIENCE. VALUE. INNOVATION.

OviTex® Reinforced Tissue Matrix

Permanent Polymer

EN

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

OviTex Reinforced Tissue Matrix with Permanent Polymer is a sterile surgical mesh. Devices are composed of ovine (sheep) derived extracellular matrix (ECM) and a monofilament polypropylene which comprises up to approximately 19 g/m² of the final product.


The polypropylene is colored with less than 0.5% ([phthalocyaninato(2-)] copper).

OviTex will incorporate into the recipient tissue with associated cellular and microvascular ingrowth. The biologic ECM components provide a scaffold for cell repopulation and to aid during soft tissue repair.


The use of surgical mesh as part of hernia repair has been shown to decrease hernia recurrence rates compared to repairs without the use of surgical mesh. OviTex Permanent is available in Core, 1S and 2S configurations (See table below).

CORE	4 layers	2 textured sides
1S	6 layers	1 smooth side and 1 textured side
2S	8 layers	2 smooth sides

The device consists of textured and/or smooth sides. The textured side of the device, indicated with blue polypropylene, provides a surface conducive for tissue ingrowth. The smooth side of the device, containing clear polypropylene, provides a surface designed to minimize tissue attachment. Image examples below

Textured side

Smooth side

OviTex is provided in various sizes to suit surgeon preference and the complexity of the soft tissue repair. The device may be trimmed to a desired shape to further accommodate an individual patient's requirements.

INDICATIONS FOR USE

OviTex is intended for use as a surgical mesh to reinforce and/or repair soft tissue where weakness exists. Indications for use include the repair of hernias and/or abdominal wall defects that require the use of reinforcing or bridging material to obtain the desired surgical outcome.

CONTRAINDICATIONS

Do not use OviTex in patients with a known sensitivity to materials of ovine (sheep) origin. Use of OviTex in this patient population may result in an allergic or immunological reaction.

RISK STATEMENTS

WARNINGS

- Device is terminally sterilized using ethylene oxide. Inspect the packaging to ensure it is intact and undamaged prior to use. Do not use OviTex, and discard prior to use if the packaging is damaged. Do not resterilize.
- Single use only. Reuse, resterilization, reprocessing and/or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device, which may result in device failure and/or patient injury. Open and unused material should be discarded.

PRECAUTIONS

- The device should not be used in patients with a known sensitivity to polypropylene.
- Do not use the product past its expiration date. The expiration date is displayed on the product labeling as the year (4 digits), month (2 digits), and day (2 digits) next to an hourglass symbol.
- Place the device in maximum contact with healthy, well-vascularized tissue to promote cell ingrowth and tissue remodeling.

POTENTIAL COMPLICATIONS

The following adverse events have been reported for surgical repair of hernias (with or without a surgical mesh): pain, infection, dysphagia, hernia recurrence, dehiscence, abscess, adhesion, bowel obstruction, bleeding, fistula, seroma, perforation, mesh migration, and mesh contraction.

MRI SAFETY INFORMATION

OviTex is MR Safe.

HOW SUPPLIED

OviTex is packaged in a Tyvek/film double pouch configuration.

LATEX INFORMATION

OviTex is not made with natural rubber latex.

STORAGE

OviTex should be stored at room temperature in a clean and dry area.

INSTRUCTIONS

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

- Inspect the packaging to ensure it is intact and undamaged.
- Using aseptic technique, peel open the outer pouch and transfer the sterile inner pouch to the sterile field.
- Open the inner pouch aseptically.
- Rehydrate the device in a sufficient volume of sterile saline or sterile Lactated Ringer's solution for a minimum of 5 minutes.
- Prepare the site using standard surgical techniques.
- Using aseptic technique, trim the device to fit the site, if necessary, providing an allowance for overlap. Position the device to achieve maximum contact between the device and surrounding tissue. To facilitate cell migration and tissue ingrowth, an overlap of 3-5 cm with healthy well-vascularized tissue is suggested.

Note: If the device is cut too small for the defect, excess tension may be placed on the suture line. This can result in recurrence of the original tissue defect or development of a defect in the adjacent tissue.

7. Using aseptic technique, transfer the device to the surgical site and suture, staple, or tack into place, avoiding excess tension.

Note: In circumstances where the OviTex 1S or 2S device contacts the viscera, place the smooth side of the product in contact with the viscera. OviTex CORE is not recommended for viscera contact.

8. Complete the surgical procedure.
9. Discard any unused portions according to institutional guidelines for medical waste.

SYMBOLS GLOSSARY

Symbols contained in ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements are indicated below.		
	5.1.1	Manufacturer
	5.1.2	European Commission Authorized Representative
	-	UK Responsible Person
	-	Swiss Authorized Representative
	-	Conformité Européene
	5.1.3	Date of manufacture
	5.1.4	Use by/Expiration date
	5.1.5	Lot number
	5.1.6	Catalog number
	5.1.8	Importer
	5.2.3	Sterilized using ethylene oxide
	5.2.6	Do not resterilize

	5.2.8	Do not use if packaging is damaged and consult instructions for use
	5.2.12	Double sterile barrier system
	5.3.4	Keep dry
	5.4.2	Do not reuse
	5.4.3	Consult instructions for use
	5.4.8	Contains biological material of animal origin
	5.7.7	Medical Device
	-	Number of units
Rx Only	-	Prescription Only. Caution: Federal (USA) law restricts this device to sale by or on the order of a physician

Report any serious incident related to OviTex to TELA Bio via S4M Europe, as well as the applicable Member State competent authority.

Users within the EU may access a Summary of Safety and clinical performance (SSCP) using Basic UDI-DI 09421900VITEXZ7 through EUDAMED or upon request (<https://ec.europa.eu/tools/eudamed>).

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THIS PAGE IS INTENTIONALLY LEFT BLANK

THIS PAGE IS INTENTIONALLY LEFT BLANK

OVITEX® CORE

REINFORCED TISSUE MATRIX

OVITEX® 1S

REINFORCED TISSUE MATRIX

OVITEX® 2S

REINFORCED TISSUE MATRIX

POLYMER REINFORCEMENT

PERMANENT

Polypropylene

Distributed by:

1 Great Valley Parkway, Suite 24
Malvern, PA 19355 USA
T (844) 835-2246
F (844) 455-8707
www.telabio.com

Aroa Biosurgery Ltd.
2 Kingsford Smith Place,
Airport Oaks, Auckland 2022,
New Zealand

S4M EUROPE

59, rue Castellion
01100 Oyonnax - France

MedEnvoy

Prinses Margrietplantsoen 33
– Suite 123 2595 AM The Hague
The Netherlands

MedEnvoy

Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

MedEnvoy

Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Easy Medical Device Ltd.
125 Deansgate
Manchester, UK M3 2LH

MedEnvoy UK Limited

85, Great Portland Street,
First Floor, London,
W1W 7LT
United Kingdom

AROA™, Aroa Biosurgery™ and AROA ECM™
are trademarks of **Aroa Biosurgery Limited**.

TELA Bio and **OviTex®** are registered trademarks
of **TELA Bio, Inc.** © TELA Bio, Inc.